

2021年12月新製品のご紹介

新製品
購入時に
使える

メルマガ会員限定 30% OFFクーポン プレゼント

毎月1日配信のCST NEWSにお得なクーポンが付いてくる！

JANUARY

DECEMBER							JANUARY						
S	M	T	W	T	F	S	S	M	T	W	T	F	S
1	2	3	4	5	6	7							
8	9	10	11	12	13	14	2	3	4	5	6	7	8
15	16	17	18	19	20	21	9	10	11	12	13	14	15
22	23	24	25	26	27	28	16	17	18	19	20	21	22
29	30	31					23	24	25	26	27	28	29

5 6 7 8 9 10 11 12 13 14 15 16 17 18

毎月1日はCST NEWSの日

12月の30% OFFクーポンコード有効期限

2021年12月1日 - 2022年1月31日受注分まで

CSTメルマガ登録はこちらから！
www.cst-science.com/2021npi



2021年12月新製品

12月のCST NEWSでお届けした30% OFFクーポンコード対象製品リストです

12月の30% OFFクーポンコード有効期間：2021年12月1日-2022年1月31日受注分まで

製品番号	製品名	サイズ	希望納入 価格 (円)	クーポン適用 特別価格 (円)	アプリケーション	種交差性
Adhesion/ECM						
63034S	COL3A1 (E8D7R) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IHC-P	H
63034T		20 µl	39,000	27,300		
74978S	MUC5B (E4Y5Z) Rabbit mAb	100 µl	65,000	45,500	IHC-P, IF-IC	H
44454S	MUC13 (E6Z1K) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IHC-P, IF-IC	H
49536S	PAI-1 (E3I5H) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IHC-P	H
49536T		20 µl	39,000	27,300		
20302S	Periostin (E5F2S) Rabbit mAb	100 µl	65,000	45,500	WB, IP, IHC-P	H
Apoptosis						
40618S	Gasdermin E Antibody	100 µl	63,000	44,100	WB	M, R
88169S	Mcl-1 (D2W9E) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	66,000	46,200	IF-IC, F	H, M, R
66026S	PARP (46D11) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	F	H, M, R, Mk
53727S	PARP (46D11) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	66,000	46,200	F	H, M, R, Mk
34366S	PARP (46D11) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	66,000	46,200	F	H, M, R, Mk
58274S	Phospho-RIP (Ser320) Antibody	100 µl	72,000	50,400	WB	H
Autophagy						
80158S	CCPG1 (E3C5G) Rabbit mAb	100 µl	63,000	44,100	WB, IP	H
12074S	Phospho-SQSTM1/p62 (Ser349) (E7M1A) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	68,000	47,600	IF-IC, F	H, M, (R)
Ca, cAMP and Lipid Signaling						
37067S	PhosphoPlus® PLCγ2 (Tyr759) Antibody Duet	1 Kit	121,000	84,700	-	-
Cell Cycle / Checkpoint Control						
99274S	Phospho-Rb (Ser807) Antibody	100 µl	72,000	50,400	WB	H
76998S	XRCC1 (E4A3V) XP® Rabbit mAb	100 µl	67,000	46,900	WB, IP, IHC-P, IF-IC	H, M, R
76998T		20 µl	39,000	27,300		
Chromatin Regulation / Nuclear Function						
50607S	FTSJ3 (E6T9J) Rabbit mAb	100 µl	63,000	44,100	WB	H, M
20123S	Tri-Methyl-Histone H3 (Lys27) (C36B11) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	68,000	47,600	F	H, M, R, Mk, (X, Z)
10584S	Tri-Methyl-Histone H3 (Lys27) (C36B11) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	68,000	47,600	F	H, M, R, Mk, (X, Z)
88197C	FastScan™ LbCpf1 (Strain ND2006) ELISA Kit	1 Kit	115,000	80,500	ELISA	-
34907S	MLL1 (E9W9Z) Rabbit mAb (Carboxy-terminal Antigen)	100 µl	63,000	44,100	WB, IP, ChIP	H, Mk
79030S	TIF1α/TRIM24 (E9T3N) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IHC-P	H, M, R, Mk
92145S	TTF-I Antibody	100 µl	63,000	44,100	WB	H, M, R, Mk
Cytoskeletal Signaling						
47610S	KIF5B (E3H4J) Rabbit mAb	100 µl	63,000	44,100	WB, IF-IC	H, M, R
62696S	KIF5B (E4A5A) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IF-IC	H
47293S	SETD3 Antibody	100 µl	63,000	44,100	WB, IP	H, Mk
37430S	ZNF185 (E6W4I) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IF-IC	H
86942S	ZO-1 (D6L1E) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	IF-IC	H, Mk
98225S	ZO-1 (D6L1E) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	66,000	46,200	IF-IC	H, Mk
Developmental Biology						
64022S	CD15/SSEA1 (MC480) Mouse mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	IF-IC, F	M
16274S	ETV5 (E5G9V) Rabbit mAb	100 µl	63,000	44,100	WB, IP	H, M
54319S	Human FGF-9 Recombinant Protein	10 µg	46,000	32,200	-	-
74356S	Mouse FGF-9 Recombinant Protein	10 µg	46,000	32,200	-	-
41560S	HIF-1α (D1S7W) XP® Rabbit mAb (Biotinylated)	100 µl	76,000	53,200	WB	H, M, Mk
70948S	Notch2 (D76A6) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	53,200	IF-IC, F	H, M, R
36135S	POU2F3 (E5N2D) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IHC-Bond, IHC-P	H
36135T		20 µl	39,000	27,300		
18113S	RUNX3/AML2 (E8D5W) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IP, IHC-P	H
18113T		20 µl	39,000	27,300		
44909S	Human TGF-α Recombinant Protein	25 µg	21,000	14,700	-	-
10179S	YAP/TAZ (D24E4) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	66,000	46,200	F	H, M, Mk

製品番号	製品名	サイズ	希望納入 価格 (円)	クーポン適用 特別価格 (円)	アプリケーション	種交差性
Glucose / Energy Metabolism						
66933S	GPAT4 Antibody	100 µl	63,000	44,100	WB	H
Immunology and Inflammation						
81123S	ASC/TMS1 (D2W8U) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	66,000	46,200	F	M
29956S	ASC/TMS1 (D2W8U) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	66,000	46,200	F	M
90776S	Bcl-3 (E2M6N) Rabbit mAb	100 µl	63,000	44,100	WB, F	H
16611S	Bcl-3 (E4D4K) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IHC-P	H, M
64137S	Phospho-Btk (Tyr551) (E5Y6N) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	68,000	47,600	F	H
16803S	Phospho-Btk (Tyr551) (E5Y6N) Rabbit mAb (PE Conjugate)	100 µl	68,000	47,600	F	H
88083S	CD3ζ Antibody	100 µl	63,000	44,100	WB	H, M
23578S	CD27 Antibody	100 µl	63,000	44,100	WB	M
26233S	CD34 (E2J1K) Rabbit mAb	100 µl	63,000	44,100	WB, IF-F, IF-IC	M
46894S	CD37 (E4K2M) Rabbit mAb	100 µl	63,000	44,100	WB, IP, F	H
55084S	CD74 (D5N3I) XP® Rabbit mAb (PE Conjugate)	100 µl	76,000	53,200	F	H
89660S	Phospho-CD79A (Tyr210) (E8E9Z) Rabbit mAb	100 µl	72,000	50,400	WB, F	H
57950S	CLEC12A/CLL-1 (E5T5V) Rabbit mAb	100 µl	63,000	44,100	WB, IP, F	H
53560S	CTLA-4 (E2V1Z) Rabbit mAb	100 µl	65,000	45,500	WB, IP, IHC-Bond, IHC-P	H, M
87223S	DPEP-1/MDP (E4O6J) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IHC-P, IF-IC	H
87223T		20 µl	39,000	27,300		
76290S	DPEP-1/MDP (E5C2V) Rabbit mAb	100 µl	63,000	44,100	WB, IF-IC	H
78401S	FcRγ (E6Y1A) Rabbit mAb	100 µl	65,000	45,500	WB, IP, IHC-Bond, IHC-P, IF-F, IF-IC	H, M
20085S	FcRn (E4D5I) Rabbit mAb	100 µl	65,000	45,500	WB, IHC-P	M, R
45850S	Fibrinogen gamma chain (E1U3Z) XP® Rabbit mAb	100 µl	70,000	49,000	WB, IHC-Bond, IHC-P	H
45850T		20 µl	39,000	27,300		
89572S	Galectin-3/LGALS3 (E7B6R) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IHC-P, IF-F, IF-IC	H, M, R
99695S	Granzyme H (E3H7W) Rabbit mAb	100 µl	63,000	44,100	WB, F	H
61606S	HPK1 (E1C3L) Rabbit mAb (Biotinylated)	100 µl	66,000	46,200	WB	H
18526S	IL-2Rγ (E2Z5W) Rabbit mAb	100 µl	63,000	44,100	WB, IP, F	M
57529S	Human IL-11 Recombinant Protein	10 µg	46,000	32,200	-	-
76561S	Mouse IL-11 Recombinant Protein	10 µg	46,000	32,200	-	-
68350S	Integrin α11 Antibody	100 µl	63,000	44,100	WB, IP	H
98344S	IRF-8 (E8X4K) Rabbit mAb	100 µl	63,000	44,100	WB, IP, IF-F, IF-IC, F, ChIP	H, M, R
51817S	ITK (E4X7M) Rabbit mAb (Biotinylated)	100 µl	66,000	46,200	WB	H
89368S	LAMP1 (E5N9Z) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	IF-F, IF-IC	M
73294S	LAMP1 (E5N9Z) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	66,000	46,200	IF-F, IF-IC	M
98215S	MSR1 (E7P6W) Rabbit mAb	100 µl	65,000	45,500	IHC-P	M
10054S	NEK7 Antibody	100 µl	63,000	44,100	WB	H, M
96874S	Phospho-NF-κB p65 (Ser529) Antibody	100 µl	72,000	50,400	WB	H
84223T	Human Reactive M1 vs M2 Macrophage IHC Antibody Sampler Kit	1 Kit	117,000	81,900	-	-
97624T	Mouse Reactive M1 vs M2 Macrophage IHC Antibody Sampler Kit	1 Kit	117,000	81,900	-	-
73582C	PathScan® Phospho-Stat1 (Tyr701) Sandwich ELISA Kit	1 Kit	114,000	79,800	ELISA	H, M
73452S	Phospho-Stat3 (Tyr705) (D3A7) XP® Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	78,000	54,600	F	H, M, R, Mk, (Hm, B, Pg, Hr)
61958S	Phospho-Stat3 (Tyr705) (D3A7) XP® Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	78,000	54,600	F	H, M, R, Mk, (Hm, B, Pg, Hr)
72337S	TNF-R2 (E8D7P) Rabbit mAb	100 µl	63,000	44,100	WB, IP, F	H
38519S	Toll-like Receptor 4 (E5D8T) Rabbit mAb	100 µl	63,000	44,100	WB	H
58125S	Tox (E5J8A) Rabbit mAb	100 µl	63,000	44,100	IF-IC, F	H
Infectious Diseases						
74698S	SARS-CoV-2 Envelope Protein Antibody	100 µl	63,000	44,100	WB	Vir
37393C	PathScan® SARS-CoV-2 Nucleocapsid Protein Sandwich ELISA Kit	1 Kit	114,000	79,800	ELISA	Vir
34340S	SARS-CoV-2 ORF3a Antibody	100 µl	63,000	44,100	WB, IP	Vir
MAP Kinase Signaling						
40737S	MEK5 (E4C4B) Rabbit mAb	100 µl	63,000	44,100	WB, IP	H, M, R
Metabolism						
58662S	AGK Antibody	100 µl	63,000	44,100	WB	H
24972S	AlfM2/FSP1 Antibody	100 µl	63,000	44,100	WB	H
31219S	COX2/MT-CO2 Antibody	100 µl	63,000	44,100	WB, IP	H, Mk
26381S	DHODH (E9X8R) Rabbit mAb	100 µl	63,000	44,100	WB	H, M, R
17412S	MT-ATP8 Antibody	100 µl	63,000	44,100	WB, IP	H
45535S	Perilipin-2 Antibody	100 µl	63,000	44,100	WB, IP	H

製品番号	製品名	サイズ	希望納入 価格 (円)	クーポン適用 特別価格 (円)	アプリケーション	種交差性
65675S	TOM70 Antibody	100 µl	63,000	44,100	WB	H, M, R
53370S	Xanthine Dehydrogenase/Oxidase Antibody	100 µl	63,000	44,100	WB	H

Motif Antibodies

46784S	PTMScan® HS Acetyl-Lysine Motif (Ac-K) Kit	1 Kit	お問い合わせ	対象外	-	-
--------	--	-------	--------	-----	---	---

Neuroscience

44864S	Phospho-CREB (Ser133) (87G3) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	68,000	47,600	F	H, M, R, (Z)
54960S	Phospho-CREB (Ser133) (87G3) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	68,000	47,600	F	H, M, R, (Z)
10767S	EAAT2 (E3P5K) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	IF-F	H, M, R
18446S	GAT1 (E7J1B) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	66,000	46,200	IF-F	H, M, R
73561S	ITM2A (E7G8A) Rabbit mAb	100 µl	63,000	44,100	WB, IF-IC	H, (M, R)
59576S	Phospho-LRRK2 (Ser935) Antibody	100 µl	72,000	50,400	WB	H, M
65915S	...	100 µl	67,000	46,900	WB, IP, IHC-P, IF-F	H, M, R
65915T	Olig2 (E6G6Q) XP® Rabbit mAb	20 µl	39,000	27,300		
32833S	Parkin (D4Z1W) Rabbit mAb	100 µl	63,000	44,100	WB	H, M, R
46421S	Phospho-PINK1 (Ser228) Antibody	100 µl	72,000	50,400	WB	H, M, (R)
56016S	Stathmin-2 Antibody	100 µl	63,000	44,100	WB, IP	H, M, R
96767S	Tau (Tau46) Mouse mAb (Biotinylated)	100 µl	66,000	46,200	WB	H, M, R, (B)
29443C	PathScan® Tau 4R Sandwich ELISA Kit	1 Kit	114,000	79,800	ELISA	H, M, R
62672S	Phospho-Tau (Thr181) (D9F4G) Rabbit mAb (Biotinylated)	100 µl	68,000	47,600	WB	H, M, R
20473S	Phospho-Tau (Thr231/Ser235) (E7G5W) Rabbit mAb	100 µl	72,000	50,400	WB, IP	H, M, R

Nuclear Receptor Signaling

20597S	Androgen Receptor (AR-V7 Specific) (E3O8L) Rabbit mAb (Biotinylated)	100 µl	66,000	46,200	WB	H
96428S	NR2F6/EAR-2 Antibody	100 µl	63,000	44,100	WB	H, M, R
74076S	PPARδ (E8O3H) Rabbit mAb	100 µl	63,000	44,100	WB, IP, ChIP	H, M, R, Mk

PI3K / Akt Signaling

77440S	Phospho-Akt (Thr308) (D25E6) XP® Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	78,000	54,600	F	H, M, R, Mk
--------	--	--------	--------	--------	---	-------------

Protein Folding and Trafficking

38700S	Transportin-1 Antibody	100 µl	63,000	44,100	WB	H, M, R, Mk
--------	------------------------	--------	--------	--------	----	-------------

Translational Control

73957S	eIF2γ Antibody	100 µl	63,000	44,100	WB	H, M, R, Mk
15719S	PTBP2 Antibody	100 µl	63,000	44,100	WB	H, M, R, Mk
57540S	TRMT112 Antibody	100 µl	63,000	44,100	WB	H, M, R, Mk

Ubiquitin and Ubiquitin-like proteins

28384C	PathScan® Phospho-Ubiquitin (Ser65) Sandwich ELISA Kit	1 Kit	114,000	79,800	ELISA	H
87525S	Phospho-VCP (Ser13) Antibody	100 µl	72,000	50,400	WB	H, (M, R)

Related Products

15764S	SignalStain® Boost IHC Detection Reagent (AP, Rat)	15 ml	68,000	47,600	IHC-P	-
26927S	SignalStain® Boost IHC Detection Reagent (AP, Goat)	15 ml	68,000	47,600	IHC-P	-
37297S	HA-Tag (C29F4) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	66,000	46,200	IF-IC, F	All
31628S	V5-Tag (D3H8Q) Rabbit mAb (Magnetic Bead Conjugate)	400 µl	66,000	46,200	IP	All
53475S	V5-Tag (E9H8O) Mouse mAb (Magnetic Bead Conjugate)	400 µl	66,000	46,200	IP	All

種交差性：Hヒト、Mマウス、Rラット、Hm-/Hmスター、Mkサル、Mkミンク、Cニワトリ、Pgブタ、Bウシ、Dmキイロショウジョウバエ、Zゼブラフィッシュ、Xアフリカツメガエル、Sc出芽酵母、Hrウマ、Qe線虫、GPモルモット、Dgイヌ、Rabウサギ、Vrウイルス、Altすべて
括弧付きの動物種は、配列が100%相同であるため反応すると推定されます。種交差性は、ウェスタンブロットングで確認しています。
アプリケーション：WBウェスタンブロットング、IP免疫沈降、IF-IC免疫蛍光細胞染色、IF-P免疫蛍光染色(バラフィン切片)、IF-F免疫蛍光染色(凍結切片)、IHC-P免疫組織化学染色(バラフィン切片)、IHC-Bo nd-Leica® Bond™使用免疫組織化学染色(バラフィン切片)、
IHC-F免疫組織化学染色(凍結切片)、F-フローサイトメトリー、ChIP-クロマチン免疫沈降、ChIP-seqクロマチン免疫沈降シーケンス、C&R-GIT&Hm、ELISA-ELISA
※種交差性とアプリケーションに関する最新情報は、ウェブをご覧ください。

キャンペーンの条件：本キャンペーンは日本国内においてのみ有効です。12月発売の新製品 (対象製品についてはwww.cst-science.com/2021npiをご参照ください) につき、定価の30%を割引いたします。有効期限は2022年1月31日です。現金または現金同等物、代替品の提供はありません。本キャンペーンは、新規の注文のみに有効であり、制限が適用されます。本キャンペーンは、既存の注文、保留中の注文、または過去の注文には適用されません。本キャンペーンは、キャンペーン用パンフレット、およびwww.cst-science.com/2021npiに掲載された12月発売の新製品の購入すべてに適用されます。その他のすべての製品/サービス (PTMScan®サービスを含む) およびカスタム注文は対象外です。割引の適用には、注文時にメールで受け取ったクーポンコードの記載が必要です。クーポンコードの記載がない注文は、通常価格が適用されます。他のプロモーションや割引との併用はできません。条件を満たさない製品、税、送料、手数料は、製品の小計の計算には含まれません。複製または譲渡された場合および法律で禁止される場合は無効となります。CSTは、いかなる時点においても事前の予告なく本キャンペーンを中止または変更する権利を留保します。

本パンフレット中の情報は2021年11月現在のもので、予告なく変更される場合がありますので予めご了承ください。
© 2021 Cell Signaling Technology, Inc. Cell Signaling Technology, CST, PathScan, FastScan, PTMScan, SignalStain and XP are trademarks of Cell Signaling Technology, Inc. All other trademarks are the property of their respective owners.

試験研究用 21-LRG-333847-FLY1-J1

■ 輸入販売元



〒101-0047
東京都千代田区内神田 1-6-10 笠原ビルディング 10 階
www.cellsignal.jp

■ 取扱店