

2023年1月新製品

CSTの高度に検証された 信頼ある新製品を使いませんか？

CST研究者が自信を持ってお届けする新製品が続々登場しています。
毎月1日は新製品をチェックしましょう！


今月のおすすめ製品

BATF3 (E3F3L) Rabbit mAb #72565

樹状細胞の発生に不可欠な転写因子

BATF3 (Basic leucine zipper transcriptional factor ATF-like 3, p21SNFT) は、塩基性ロイシンジッパーファミリーのメンバーであり、Jun タンパク質と二量体を形成し、AP-1 (Activator protein-1) の活性を抑制します。BATF3は、従来型の1型樹状細胞 (conventional DC1:cDC1) の発生に不可欠です。cDC1は、ナイーブCD8陽性T細胞に対してMHCクラスIを介して抗原をクロスプレゼンテーションすることにより、活性化細胞傷害性CD8陽性T細胞への分化誘導を行います。

キーワード: 樹状細胞、発生、免疫



IHC-P


IHC analysis of paraffin-embedded human Hodgkin lymphoma using #72565.

Cathepsin K (E7U5N) Rabbit mAb #57056

骨代謝に重要な酵素

Cathepsin K (CTSK) は、共通のパパイン様折り畳み構造と、保存された活性部位であるCys-Asn-Hisの3残基をもつカテプシンファミリーに属する、システインプロテアーゼです。Cathepsin Kはまず、不活性なCathepsin K前駆体として合成され、続くプロセッシングにより成熟酵素となることが知られており、その最も研究されている機能は「骨吸収の媒介」です。多くの製薬会社が、抗骨粗鬆症薬として、Cathepsin Kの選択的阻害剤の開発に取り組んでいます。また、Cathepsin Kは最初に破骨細胞で発見されましたが、Cathepsin Kの発現は、前立腺がんや乳がん、骨がん、肺がん、腎臓がんなどの様々ながんにも関与していることが分かっています。


キーワード: 骨代謝、骨粗鬆症、がん




IHC-P

IHC analysis of paraffin-embedded 4T1 syngeneic tumor using #57056.


詳細はこちらをご覧ください:
www.cst-science.com/2023npi



GPVI (E4Y9J) Rabbit mAb #45649






Phospho-p53 (Ser15) (16G8) Mouse mAb (BSA and Azide Free) #56979



TAZ (E9J5A) XP® Rabbit mAb (BSA and Azide Free) #49756

製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性	ピックアップ 製品	キャリアフリー 抗体
Adhesion/ECM						
93790T EpCAM (E6V8Y) XP® Rabbit mAb (Mouse Preferred)	20 µl	39,000	WB, IP, IHC-LB, IHC-P	H, M, R		
38313T GPNMB (E4D7P) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P	H		
45649S GPVI (E4Y9J) Rabbit mAb	100 µl	68,000	WB	H	●	
30075S MMP-1 (E8G4F) Rabbit mAb	100 µl	68,000	WB, IF-IC	H		
16564T MUC1-C (D5K9I) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P, IF-IC	H		
56572S MXRA8 (E7K5Q) Rabbit mAb	100 µl	68,000	WB	H, Mk		
25426SF Olfm4 (D6Y5A) XP® Rabbit mAb (Mouse Specific) (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F	M		●
Apoptosis						
97774S Cleaved Caspase-3 (Asp175) (D3E9) Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	80,000	FC-FP	H, (R, B, Dg, Pg)		
24844S Cleaved Caspase-7 (Asp198) (D6H1) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	80,000	FC-FP	H, M, R		
45330S Cleaved Caspase-7 (Asp198) (D6H1) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	80,000	FC-FP	H, M, R		
40446SF RIP (D94C12) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC, FC-FP	H, M, R, Hm, Mk		●
96323SF Phospho-RIP (Ser166) (D8I3A) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IF-IC, FC-FP	H		●
Autophagy						
58139SF LC3A/B (D3U4C) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC, FC-FP	H, M, R, (X, B, Dg, Pg)		●
Ca, cAMP and Lipid Signaling						
58473S Calreticulin (D3E6) XP® Rabbit mAb (PE-Cy7® Conjugate)	100 µl	81,000	FC-FP	H, M, R		
Cell Cycle / Checkpoint Control						
60721SF 53BP1 (E7N5D) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H		●
91590T Aurora A (D3V7T) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P	H		
65173SF Cyclin B1 (D5C10) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC, FC-FP	H, R		●
54011SF MCM2 (D7G11) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H, M, R, Mk		●
56979SF Phospho-p53 (Ser15) (16G8) Mouse mAb (BSA and Azide Free)	100 µg	151,000	WB, IF-IC, FC-FP	H		●
39692T p63 (D9L7L) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P	H		
53692SF Wee1 (D10D2) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, FC-FP	H, Mk		●
Chromatin Regulation / Nuclear Function						
87855SF Cas9 (<i>S. aureus</i>) (E4G3U) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC, FC-FP	All		●
52177SF Cas9 (<i>S. aureus</i>) (6H4) Mouse mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC	All		●
28754SF DNMT1 (D63A6) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, M, R, Mk, (Hm, B, Dg, Hr, GP)		●
17705S GATAD2A (E7B5B) Rabbit mAb	100 µl	68,000	WB, IP, ChIP	H, M, R, Mk		
57796S GLTSCR1L (E4D9W) Mouse mAb	100 µl	68,000	WB, IP	H, M, R, Mk		
59141S Phospho-Histone H2A.X (Ser139) (20E3) Rabbit mAb (Pacific Blue™ Conjugate)	100 µl	80,000	FC-FP	H, M, R, Mk		
68829S Histone H3 (D1H2) XP® Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	81,000	FC-FP	H, M, R, Mk, (Hm, C, Dm, X, Z, B)		
43185S Phospho-Histone H3 (Ser10) (D2C8) XP® Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	85,000	FC-FP	H, M, R, Mk, Z		
85470SF HOXB13 (D7N8O) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, Mk		●






Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) (197G2) Rabbit mAb (BSA and Azide Free) #9524

CPS1/HepPar-1(OCH1E5)MouseAb(Clone previously known as Hepatocyte Specific Antigen) #96422

Xanthine Dehydrogenase/Oxidase (E2U8Y) Rabbit mAb #35792

製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性	ピックアップ 製品	キャリアフリー 抗体
46160S MBNL1 (E8Q7B) Rabbit mAb	100 µl	68,000	WB	H, M, R, Mk		
6891T Menin (D45B1) XP® Rabbit mAb	20 µl	39,000	WB, IF-IC	H, M, R, Mk, (B, Pg, Hr)		
Cytoskeletal Signaling						
93473SF β-Actin (13E5) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC, FC-FP	H, M, R, Mk, B, Pg, (Hm, C, Dg, Hr, Rab)		●
94369SF Keratin 20 (D9Z1Z) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H		●
90307S MYO1D (E6T6Q) Rabbit mAb	100 µl	68,000	WB	H, M, R		
Developmental Biology						
37995SF CD105/Endoglin (3A9) Mouse mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H		●
64326T CD133 (D2V8Q) XP® Rabbit mAb	20 µl	39,000	WB, IHC-LB, IHC-P, IF-IC	H		
47061T CELSR2 (D2M9H) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P	H		
15171T FAM3C (D1S2D) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P, IF-IC, FC-FP	H		
3695T LIN28A (D84C11) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P, IF-IC, FC-FP	H, (R, Mk)		
3608T Notch1 (D1E11) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P, ChIP	H, M, R		
60433T Pax6 (D3A9V) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P, IF-IC	H, (M, R)		
88386SF Sox10 (E6B61) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P	H		●
49756SF TAZ (E9J5A) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P	H, M, Mk		●
Glucose / Energy Metabolism						
70720T ASS1 (D4O4B) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P, IF-IC	H, M, R		
Immunology and Inflammation						
86772S Rapid-Act T Cell Activation Kit (Mouse, Anti-CD3/CD28)	1 Kit	76,000	FUNC	M		
97279SF Aiolos (D1C1E) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, FC-FP	H, M		●
17507T ASC/TMS1 (D2W8U) Rabbit mAb (Mouse Specific) (Alexa Fluor® 488 Conjugate)	20 µl	39,000	IF-IC, FC-FP	M		
72565S BATF3 (E3F3L) Rabbit mAb	100 µl	70,000	WB, IP, IHC-P	H	●	
60066S BST2 (E5D1Z) Rabbit mAb	100 µl	68,000	WB	M		
76536S CARD14 (E3F8E) Rabbit mAb	100 µl	68,000	WB	H		
57056S Cathepsin K (E7U5N) Rabbit mAb	100 µl	70,000	WB, IHC-P	M	●	
13333T CD79A (D1X5C) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P	H		
69659SF Evi-1 (C50E12) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC, FC-FP	H		●
87985T Galectin-3/LGALS3 (D4I2R) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P, IF-IC	H		
3892T HS1 (D5A9) XP® Rabbit mAb (Rodent Specific)	20 µl	39,000	WB, IP, IHC-P, IF-F, IF-IC, FC-FP	M, R		
89799S Granzyme C Antibody	100 µl	66,000	WB	M		
16433S Granzyme F Antibody	100 µl	66,000	WB	M		
30979SF Ikaros (D6N9Y) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC, FC-FP	H, M		●
36214SF Phospho-IKKα/β (Ser176/180) (16A6) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P, FC-FP	H, M, R, Mk, (B)		●
10327T Phospho-IRF-3 (Ser396) (D6O1M) Rabbit mAb (Alexa Fluor® 647 Conjugate)	20 µl	39,000	FC-FP	H, M, R		
73236SF MLANA/MART-1 (E9Q4O) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H		●
49422S Myeloperoxidase (E1E7I) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	81,000	IF-IC, FC-FP	H		
17646S Myeloperoxidase (E1E7I) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IF-IC, FC-FP	H		




c-Jun (60A8) Rabbit mAb (BSA and Azide Free)
#26325


PTMScan® Nitro-Tyrosine Motif [Nitro-Y] Kit #44462

Late-Onset Alzheimer's Disease Risk Gene (Mouse Model) Antibody Sampler Kit #97888


製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性	ピックアップ 製品	キャリアフリー 抗体
19360S PGLYRP-1/PGRP-S (E3B6M) Rabbit mAb	100 µl	68,000	FC-FP	H		
16318T Phospho-SLP-76 (Ser376) (D7S1K) Rabbit mAb (PE Conjugate)	20 µl	39,000	FC-FP	H		
90145S Phospho-Stat1 (Tyr701) (58D6) Rabbit mAb (PE-Cy7® Conjugate)	100 µl	80,000	FC-FP	H, M		
74234SF Stat3 (D1B2J) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H, M, R		●
51865SF Phospho-STING (Ser365) (DIC4T) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	IF-IC	M		●
20399S Phospho-TBK1/NAK (Ser172) (D52C2) XP® Rabbit mAb (Alexa Fluor® 532 Conjugate)	100 µl	85,000	FC-FP	H, M, (R, Mk, X, B, Dg)		
73382T Phospho-Zap-70 (Tyr319)/Syk (Tyr352) (65E4) Rabbit mAb (Alexa Fluor® 488 Conjugate)	20 µl	39,000	FC-FP	H, M, (R, Hm, C, B, Dg, Pg, Hr)		
Infectious Diseases						
95561S SARS-CoV-2 Nsp9 (E3E3X) Rabbit mAb	100 µl	68,000	WB	Vir		
MAP Kinase Signaling						
26325SF c-Jun (60A8) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, M, R, Mk		●
9524SF Phospho-p44/42 MAPK (Erk1/2) (Thr202/Tyr204) (197G2) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IF-IC, FC-FP	H, M, R, Mk, Mi, Dm, Z, Pg		●
10242SF Src (36D10) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC, FC-FP	H, M, (C)		●
Metabolism						
47213T 4F2hc/CD98 (D3F9D) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P, IF-IC, FC-FP, FC-L	H		
94274T 4F2hc/CD98 (D5U4G) XP® Rabbit mAb	20 µl	39,000	WB, IHC-P, IF-IC	H		
40487S ACBP/DBI (E4V8V) XP® Rabbit mAb	100 µl	75,000	WB, IHC-P	H, M, R		
40487T	20 µl	39,000				
45365S ACBP/DBI (E6K8G) XP® Rabbit mAb	100 µl	75,000	WB, IHC-P	H		
45365T	20 µl	39,000				
69974SF ASCT2 (D7C12) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, Mk		●
96422S CPS1/Hep Par-1 (OCH1E5) Mouse mAb (Clone previously known as Hepatocyte Specific Antigen)	100 µl	70,000	IHC-P	H	●	
94004T CYP17A1 (E6A7G) XP® Rabbit mAb	20 µl	39,000	WB, IP, IF-F	H, M, R		
92261SF Legumain (D6S4H) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, M, R		●
35939T Malic Enzyme 2 (E1N3F) XP® Rabbit mAb	20 µl	39,000	WB, IF-IC, FC-FP	H, Mk		
68533S PAPSS1 (E6U5C) Rabbit mAb	100 µl	68,000	WB	H, M, R, Mk		
95251SF PKM2 (D78A4) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, M, R, Mk		●
35792S Xanthine Dehydrogenase/Oxidase (E2U8Y) Rabbit mAb	100 µl	68,000	WB, IF-IC	H	●	
Motif Antibodies						
34200S PTMScan® Control Peptides O-GlcNAc	100 µl	48,000	-	-		
44462S PTMScan® Nitro-Tyrosine Motif [Nitro-Y] Kit	1 Kit	お問い合わせ	-	-	●	
Neuroscience						
59678T AQP4 (D1F8E) XP® Rabbit mAb	20 µl	39,000	WB, IP, IHC-P, IF-F	H, M, R		
80063SF Glutamate Dehydrogenase 1/2 (D9F7P) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H, M, R, Mk		●
84037S Iba1/AIF-1 (E4O4W) XP® Rabbit mAb (Pacific Blue™ Conjugate)	100 µl	81,000	FC-FP	H, M, R		
97888T Late-Onset Alzheimer's Disease Risk Gene (Mouse Model) Antibody Sampler Kit	1 Kit	148,000	-	M	●	



FoxO1 (C29H4) Rabbit mAb (BSA and Azide Free)
#29336



Human Reactive Exosome Marker Antibody Sampler
Kit #80610



3X DYKDDDDK Tag (E7C5T) Rabbit mAb #87537

製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性	ピッキングアップ 製品	キャリアフリー 抗体
35380SF TMEM119 (E3E1O) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	IF-F	M, (H)		●
55739S TREM2 (E4J7A) Rabbit mAb	100 µl	68,000	WB	H, M		
PI3K / Akt Signaling						
29336SF FoxO1 (C29H4) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, M, R, Mk		●
Protein Folding and Trafficking						
64293SF Phospho-HSP27 (Ser82) (D1H2F6) XP® Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P, IF-IC, FC-FP	H, M, (R, Hm, B, Dg, Hr)		●
80610T Human Reactive Exosome Marker Antibody Sampler Kit	1 Kit	138,000	-	H	●	
79842S ZDHHC5 (E7T4N) Rabbit mAb	100 µl	68,000	WB, IP	H		
RTK						
2085T EGF Receptor (E746-A750del Specific) (D6B6) XP® Rabbit mAb	20 µl	49,000	WB, IP, IHC-P, IF-P, IF-IC, FC-FP	H		
76443S EGF Receptor (D38B1) XP® Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	81,000	FC-FP	H, M, Mk		
47737SF Phospho-EGF Receptor (Tyr1173) (53A5) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P	H, M, R		●
94065SF FGF Receptor 3 (C51F2) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H		●
23991S Phospho-FLT3 (Tyr589/591) (E4D6Y) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	80,000	FC-FP	H		
Translational Control						
99623S ADAR1 (E6X9R) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IF-IC, FC-FP	H, Mk		
59278S DDX21 (E9O2Q) Rabbit mAb	100 µl	68,000	WB, IF-IC	H, M, R, Mk		
62788S Phospho-S6 Ribosomal Protein (Ser235/236) (D57.2.2E) XP® Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	85,000	FC-FP	H, M, R, Mk, Mi, Sc, (C, Pg)		
88441S Phospho-S6 Ribosomal Protein (Ser240/244) (D68F8) XP® Rabbit mAb (Alexa Fluor® 700 Conjugate)	100 µl	85,000	FC-FP	H, M, R, Mk, (C, Pg)		
18791S TRMT6 (E7L3M) Rabbit mAb	100 µl	68,000	WB, IP	H, R, Mk		
30894S ZFP36L1 (E6L6S) Rabbit mAb	100 µl	68,000	WB	H, M, Mk		
Ubiquitin and Ubiquitin-like proteins						
5591T SENP3 (D20A1O) XP® Rabbit mAb	20 µl	39,000	WB, IP, IF-IC	H, M, R, Mk		
83816S Phospho-Ubiquitin (Ser65) (E5T1W) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	80,000	FC-FP	H, (M, R)		
Related Products						
87537S 3X DYKDDDDK Tag (E7C5T) Rabbit mAb	100 µl	68,000	WB, IP, IF-IC	All	●	
20818S HA-Tag (C29F4) Rabbit mAb (Alexa Fluor® 750 Conjugate)	100 µl	76,000	FC-FP	All		

CST最新情報




CUT&RUNで実験を迅速化

CUT&RUNは比較的新しいアプリケーションのため、検証済みの市販の抗体は限られており、ChIPまたはChIP-seq検証済み抗体の使用を推奨する企業もあります。しかし、これらの抗体はスタート地点としては良いのですが、必ずしも成功を保証するものではありません。CUT&RUN実験で、ChIPやChIP-seq用の抗体が必ずしも機能しない理由についてご紹介します。

ブログはこちら：

www.cst-science.com/23npi-jan-info1




汎用性の高いCAR検出試薬の導入：抗CARリンカー抗体

キメラ抗原受容体 (CAR) T細胞療法とは、自身の免疫系を活用して抗腫瘍効果を誘導する方法です。CSTは、フローサイトメトリーパネルへ組み込みが可能で、抗原の特異性に関係なく、CARの同定および発現を確認できる汎用性の高い市場初の抗体を開発しました。

ブログはこちら：

www.cst-science.com/23npi-jan-info2



Solutions for Neurodegeneration Drug Discovery

CSTの幅広い抗体製品ラインアップや主要な指標、ready-to-use ELISA、細胞解析用キットをまとめたSolutions for Advancing Neurodegeneration Therapeutics Discoveryをダウンロードして、効果的な新規治療法の開発にお役立てください。

詳細はこちら：

www.cst-science.com/23npi-jan-info3



CSTテクニカルセミナーシリーズ

無料オンラインセミナーを毎月実施しています。弊社の技術者が培った、実験をスムーズに進めるためのコツをお伝えします。どなたでもお気軽にご参加ください。

詳細はこちら：

www.cst-science.com/23npi-jan-info4



種差特性：H-ヒト、M-マウス、R-ラット、Hm-ハムスター、Mk-サル、Mi-ミンク、C-ニワトリ、Pg-ブタ、B-ウシ、Dm-キイロショウジョウバエ、Z-ゼブラフィッシュ、X-アカツメガエル、Sc-出芽酵母、Hr-ウマ、Ce-線虫、GP-モルモット、Dg-イヌ、Rab-ウサギ、Vir-ウイルス、All-すべて
括弧付きの動物種は、配列が100%相同であるため反応すると推定されます。種差特性は、ウェスタンブロッティングで確認しています。
アプリケーション：WB-ウェスタンブロッティング、IP-免疫沈降、IF-IC-免疫蛍光細胞染色、IF-P-免疫蛍光染色 (パラフィン切片)、IF-F-免疫蛍光染色 (凍結切片)、IHC-P-免疫組織化学染色 (パラフィン切片)、IHC-LB-IHC for Leica BOND platform 使用免疫組織化学染色 (パラフィン切片)、IHC-F-免疫組織化学染色 (凍結切片)、FC-FP-フローサイトメトリー (固定、透過化)、FC-L-フローサイトメトリー (生細胞)、ChIP-クロマチン免疫沈降、ChIP-seq-クロマチン免疫沈降シーケンス、C&R-CUT&RUN、ELISA-ELISA、R Dot Blot-RNAドットプロット、FUNC-Functional

本パンフレット中の情報は2022年12月現在のもので、予告なく変更される場合がありますので予めご了承ください。

© 2022 Cell Signaling Technology, Inc. Cell Signaling Technology, CST, PTMScan, and XP are trademarks or registered trademarks of Cell Signaling Technology, Inc. Alexa Fluor® is a registered trademark of Life Technologies Corporation. Pacific Blue™ is a trademark of Life Technologies Corporation. All other trademarks are the property of their respective owners. Visit cellsignal.com/trademarks for more information.

輸入販売元



〒101-0047

東京都千代田区内神田1-6-10 笠原ビルディング10階

www.cellsignal.jp

取扱店