

2023年7月新製品

CSTの高度に検証された信頼ある新製品を使いませんか？



CST研究者が自信を持ってお届けする新製品が続々登場しています。
毎月1日は新製品をチェックしましょう！

今月のおすすめ製品


Clusterin (E2O1D) Rabbit mAb #23159

分子シャペロン

Clusterinは、ストレス誘導性タンパク質と相互作用し、安定化させて蓄積を防ぐ、分泌シャペロンタンパク質として機能します。アルツハイマー病においては、アミロイドβ (1-40) ペプチドと長期間にわたる安定した複合体を形成して隔離し、アミロイド線維の形成を阻止することにより、保護的な役割を果たしていることが研究から分かっています。

また、分泌タンパク質としてだけでなく、いくつかのアイソフォームが、細胞内の核やミトコンドリア、細胞質、小胞体に局在しています。これらの複数のアイソフォームの細胞内分布が、Clusterinの機能の多様性につながっています。さらに、Clusterinは、膜のリサイクリングや細胞接着、細胞増殖、アポトーシス、腫瘍の生存に関与していることが報告されています。

キーワード：分子シャペロン、アルツハイマー病



IF-F

Confocal IF analysis of fixed frozen mouse hippocampus labeled with #23159 (left, green). Free secondary binding sites were then blocked with Rabbit (DA1E) mAb IgG XP® Isotype Control #3900 prior to co-labeling with GFAP (GA5) Mouse mAb (Alexa Fluor® 647 Conjugate) #3657 (right, red), Iba1/AIF-1 (E4O4W) XP® Rabbit mAb (Alexa Fluor® 555 Conjugate) #36618 (right, cyan pseudocolor), and ProLong Gold Antifade Reagent with DAPI #8961 (right, blue).


BCL2L13 (E6G1Q) Rabbit mAb #61974

アポトーシスとミトファジーを制御

BCL2L13は、ミトコンドリアにおいてアポトーシスを制御するBcl-2ファミリータンパク質の非定型的なメンバーです。BCL2L13は、細胞内でアポトーシスを制御している可能性があります。また、BCL2L13はアデニンヌクレオチドトランスロカーゼ (ANT) と相互作用し、ANTのADP/ATP依存性の交換輸送活性を抑制します。他の研究では、BCL2L13が膠芽細胞腫におけるアポトーシスを阻害し、腫瘍の成長を促進することが報告されています。これらの研究から、BCL2L13はアポトーシス促進性のCerS2とCerS6に結合し、CerS2/6のホモ複合体やヘテロ複合体の形成および活性を阻害することが分かりました。

そして、BCL2L13はミトコンドリアの断片化を調節し、ミトファジーを介して傷害を受けたミトコンドリアの除去を促進するためにLC3に結合することが分かっています。

キーワード：アポトーシス、ミトファジー、細胞死



WB


WB analysis of extracts from various cell lines using #61974 (upper) or GAPDH (D16H11) XP® Rabbit mAb #5174 (lower). Negative expression of BCL2L13 protein in OVSARO cells is consistent with the predicted expression pattern.

For Research Use Only. Not for Use in Diagnostic Procedures.


詳細はこちらをご覧ください: cst-science.com/2023npi




Cell Signaling
TECHNOLOGY®



Apoptosis/Necroptosis Antibody Sampler Kit II
#73966




Chk1 (E9M4D) Rabbit mAb #37010




Mismatch DNA Repair (MMR) Antibody Sampler Kit
#89884


製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性
Adhesion/ECM				
61408S DPP4/CD26 (E2V5I) Rabbit mAb	100 µl	72,000	WB	M, R
53884S LAMC2 (E9F7M) Rabbit mAb	100 µl	72,000	WB	H
47449S MMP-9 (E3D1W) Rabbit mAb	100 µl	72,000	WB	H
67950C PathScan® RP Total FAK Sandwich ELISA Kit	1 Kit	126,000	ELISA	H, M
86460S TUFT1 (E9M2U) Rabbit mAb	100 µl	72,000	WB	H
Apoptosis				
73966T Apoptosis/Necroptosis Antibody Sampler Kit II	1 Kit	138,000	-	-
61974S BCL2L13 (E6G1Q) Rabbit mAb	100 µl	72,000	WB, IP	H, M, R
16046S Cleaved Caspase-6 (Asp179) (E8L9T) Rabbit mAb	100 µl	79,000	WB	H
26776T Gasdermin Family Antibody Sampler Kit	1 Kit	117,000	-	-
Ca, cAMP and Lipid Signaling				
87716SF Calretinin (E7R6O) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC	H, M, R
Cell Cycle / Checkpoint Control				
85785S Anillin (E9L6E) Rabbit mAb	100 µl	72,000	WB	H
37010S Chk1 (E9M4D) Rabbit mAb	100 µl	72,000	WB	H, M, R
58268SF DNA-PKcs (E6U3A) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H
63273S Ki-67 (D2H10) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	IHC-P	H
88338S Ki-67 (D2H10) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	IHC-P	H
89884T Mismatch DNA Repair (MMR) Antibody Sampler Kit	1 Kit	107,000	-	-
37698SF MSH2 (D24B5) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H, M
86714S XRCC2 (E7M8Y) Rabbit mAb	100 µl	72,000	WB, IP	H
Chromatin Regulation / Nuclear Function				
62935S BAF Antibody	100 µl	66,000	WB	H, M, R, Mk
93589T BAF Complex IHC Antibody Sampler Kit	1 Kit	154,000	-	-
68280S DGCR8 (E9O3L) Rabbit mAb	100 µl	72,000	WB, IP, ChIP	H, M, R, Mk
40653S PHF6 (E8F2L) Rabbit mAb	100 µl	72,000	WB, IP, ChIP	H, M, R, Mk
37315S RBM9 (E7W6U) Rabbit mAb	100 µl	72,000	WB, IP	H, Mk
27106C PathScan® RP BRD9 Sandwich ELISA Kit	1 Kit	126,000	ELISA	H
20055S SET1B (E7D1G) Rabbit mAb	100 µl	72,000	WB, IP, ChIP	H
Cytoskeletal Signaling				
44873S Filamin A (E8P8M) Rabbit mAb	100 µl	75,000	WB, IHC-P, IF-F	H, M, R, Mk
Developmental Biology				
33699SF HES1 (D6P2U) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, M, R, Mk
32253SF NKX3.1 (D6D2Z) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H
61048S TACE (E8R8M) Rabbit mAb	100 µl	72,000	WB, IP	H, M, R
Immunology and Inflammation				
63178S CD3ε (D7A6E™) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IHC-P	H
47434S CD3ζ (E8R1Q) Rabbit mAb	100 µl	72,000	WB, IP	H
37217S CD6 (E9Q1V) Rabbit mAb	100 µl	72,000	WB, IHC-P	M
77909S CD8α (D8A8Y) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	IHC-P	H, Mk



CD6 (E9Q1V) Rabbit mAb #37217






CD47 (E2V9V) Rabbit mAb #36096



G4S Linker (E7O2V) Rabbit mAb (PE Conjugate) #38907

製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性
Immunology and Inflammation (cont.)				
48893S CD11b/ITGAM (E3J2F) Rabbit mAb	100 µl	75,000	WB, IHC-P, IF-IC	H
37940S CD11b/ITGAM (D6XIN) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	IHC-P	H
19246S CD11b/ITGAM (D6XIN) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	IHC-P	H
79750S CD11b/ITGAM (D6XIN) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	IHC-P	H, Mk
95103S CD11c (D3V1E) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	81,000	IHC-P	H
77882S CD11c (D3V1E) XP® Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	81,000	IHC-P	H
42756S CD11c (D3V1E) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IHC-P	H
86076S CD45 (Intracellular Domain) (D9M8I) XP® Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	81,000	IHC-P	H, Mk
19744S CD45 (D9M8I) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IHC-P	H
36096S CD47 (E2V9V) Rabbit mAb	100 µl	75,000	WB, IHC-P	M
58225S CD62L/L-Selectin (E7X9I) Rabbit mAb	100 µl	75,000	WB, IHC-P	H
76220S CD206/MRC1 (E6T5J) XP® Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	81,000	IHC-P, IF-F	H, M, R
43920S CD206/MRC1 (E6T5J) XP® Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	81,000	IHC-P, IF-F	H, M, R
98492SF FcRγ (E6Y1A) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P, IF-F, IF-IC	H, M
33461SF FcRn (E9J8W) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P	H
38915SF Fibrinogen gamma chain (E1U3Z) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P	H
71645L G4S Linker (E7O2V) Rabbit mAb	300 µl	204,000	FC-L	All
71645T	20 µl	39,000		
50515L G4S Linker (E7O2V) Rabbit mAb (Alexa Fluor® 488 Conjugate)	300 µl	228,000	FC-L	All
69782L G4S Linker (E7O2V) Rabbit mAb (Alexa Fluor® 647 Conjugate)	300 µl	228,000	FC-L	All
38907L G4S Linker (E7O2V) Rabbit mAb (PE Conjugate)	300 µl	228,000	FC-L	All
88084SF Galectin-1/LGALS1 (8A12) Mouse mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H
40019SF GNLY (E2T3D) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC, FC-FP	H
29268S Granzyme B (D6E9W) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	IHC-P	H, M
71535S Granzyme B (D6E9W) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	IHC-P	H, M, Mk
16483L Ibrutinib	10 mg	39,000	-	-
81361S IL-6Ra/CD126 (E5B3P) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	FC-L	H
98797S IL-6Ra/CD126 (E5B3P) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	FC-L	H
67459S Lysozyme C (E9W6E) Rabbit mAb	100 µl	75,000	WB, IHC-P	H
44238S MHC Class I (EMR8-5) Mouse mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	FC-FP	H
57549S MHC Class I (EMR8-5) Mouse mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	FC-FP	H
69717S MHC Class I (EMR8-5) Mouse mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	FC-FP	H
36845S OASL (E7W1R) Rabbit mAb	100 µl	72,000	WB	H
76700S SAMHD1 (E8V1G) Rabbit mAb	100 µl	72,000	WB, IP, FC-FP	H
20407S SerpinB1 (E4W7B) Rabbit mAb	100 µl	75,000	WB, IHC-P	H
24201S SIT Antibody	100 µl	66,000	WB	H
40284SF Phospho-Stat5 (Tyr694) (C71E5) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P, IF-IC, FC-FP	H, M, (R, B)
42383S TAPBP (E5P1F) Rabbit mAb	100 µl	72,000	WB	H
84382S TMPRSS4 (E5F4W) Rabbit mAb	100 µl	72,000	WB	H
51468S TRAF1 (E3L8A) Rabbit mAb	100 µl	72,000	WB, IP	H






Lysozyme C (E9W6E) Rabbit mAb #67459


Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 488 Conjugate) #55809

c-Fos (E217R) XP® Rabbit mAb (BSA and Azide Free) #61396


製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性
Immunology and Inflammation <i>(cont.)</i>				
57710L Whitlow/218 Linker (E3U7Q) Rabbit mAb	300 µl	204,000	FC-L	All
57710T	20 µl	39,000		
55809L Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 488 Conjugate)	300 µl	228,000	FC-L	All
69310L Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 647 Conjugate)	300 µl	228,000	FC-L	All
62405L Whitlow/218 Linker (E3U7Q) Rabbit mAb (PE Conjugate)	300 µl	228,000	FC-L	All
MAP Kinase Signaling				
47698SF β-Arrestin 1 (D7Z3W) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, M, R
60407SF B-Raf (E3T5C) Mouse mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, M, R, Mk
61396SF c-Fos (E217R) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC	H, M, R
45307SF MEK1 (D2R1O) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H, M, R, Mk, (B, Dg, Pg, GP)
64117S Phospho-SHP-2 (Tyr542) (E8D6V) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	80,000	FC-FP	H, M, R
Metabolism				
52339S GSTM1 Antibody	100 µl	66,000	WB	M, R
48159S FBP2/FBPase 2 (E9M5T) Rabbit mAb	100 µl	72,000	WB	H
84478S MAT2A (E5I7C) Rabbit mAb	100 µl	72,000	WB	H, M, R
60406S MGEA5/OGA (E9C5U) Rabbit mAb	100 µl	72,000	WB	H, M
31490SF NaPi2b/SLC34A2 (D6W2G) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-IC	H
57117SF Perilipin-1 (D1D8) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC	H, M
38244S SLC25A40 (E9C7Y) Rabbit mAb	100 µl	72,000	WB	H
Neuroscience				
10440SF β-Amyloid (D12B2) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P	H
16696S BDNF (E6Z7C) Rabbit mAb	100 µl	72,000	WB, IP	H, M
23159S Clusterin (E2O1D) Rabbit mAb	100 µl	72,000	WB, IF-F	M, R
56914SF EAAT1 (D44E2) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-F	H, M, R
83520SF EGR1 (I5F7) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-F, IF-IC	H, M, R, (B)
29642S FOXG1 (E2W8P) Rabbit mAb	100 µl	75,000	WB, IHC-P, IF-F	H, M, (R)
93972SF NeuN (E4M5P) Mouse mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-F	H, M, R
40459SF Neurofilament-H (E7Z7G) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-LB, IHC-P, IF-F, IF-IC	H, M, R
28610S NGF (E2S9D) Rabbit mAb	100 µl	72,000	WB	H, M
52675SF p35/25 (C64B10) Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F	H, M, R
44712C PathScan® Tau (GT-38) Sandwich ELISA Kit	1 Kit	126,000	ELISA	H
Nuclear Receptor Signaling				
24847S NR2F6/EAR-2 (E7Q7G) Rabbit mAb	100 µl	72,000	WB, IP	H, M
Protein Folding and Trafficking				
44989SF Caveolin-1 (D46G3) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IHC-P, IF-F, IF-IC, FC-FP	H, M, R, Hm, Mk, B, Dg
RTK				
28110SF EGF Receptor (L858R Mutant Specific) (43B2) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P, IF-P, IF-IC, FC-FP	H
28685SF FGF Receptor 4 (D3B12) XP® Rabbit mAb (BSA and Azide Free)	100 µg	134,000	WB, IF-IC	H
91415SF Phospho-HER2/ErbB2 (Tyr1221/1222) (6B12) Rabbit mAb (BSA and Azide Free)	100 µg	151,000	WB, IHC-P	H, M



SLC25A40 (E9C7Y) Rabbit mAb #38244



PathScan® Tau (GT-38) Sandwich ELISA Kit #44712



Stress Granule Marker Antibody Sampler Kit #29530

製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性
Translational Control				
39489S eIF3B Antibody	100 µl	66,000	WB	H, M, R, Mk
29530T Stress Granule Marker Antibody Sampler Kit	1 Kit	107,000	-	-
Ubiquitin and Ubiquitin-like proteins				
88503S Phospho-Ubiquitin (Ser65) (E2J6T) Rabbit mAb (Biotinylated)	100 µl	80,000	WB	H, (M, R)


市場初のCAR-T研究用CAR検出試薬：抗CARリンカーモノクローナル抗体

CST®抗CARリンカー抗体は2023年CiteAb Innovation Awardを受賞！

市場初の抗CARリンカー抗体は、抗原の特異性に関係なく、CAR-Tの同定および発現を確認できる汎用性の高い検出試薬です。CSTは、複数の異なる標識を用いた標識抗体も提供しています。また、CSTのカスタム標識サービスを利用して、カスタム標識することも可能です。全製品の一覧は以下をご覧ください。



ブログ記事はこちら：cst-science.com/23npi-jul-info5



製品名	サイズ	希望納入 価格 (円)	アプリケーション	種交差性
Immunology and Inflammation				
18839S CAR-T Cell (G4S Linker) Transduction Efficiency Flow Cytometry Panel	1 Kit	299,000	FC-L	H
35139S CAR-T Cell (Whitlow/218 Linker) Transduction Efficiency Flow Cytometry Panel	1 Kit	299,000	FC-L	H
71645S G4S Linker (E702V) Rabbit mAb	100 µl	68,000	FC-L	All
63670S G4S Linker (E702V) Rabbit mAb (BSA and Azide Free)	100 µl	134,000	FC-L	All
50515S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	FC-L	All
90841S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 532 Conjugate)	100 µl	76,000	FC-L	All
18862S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	FC-L	All
39614S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	76,000	FC-L	All
69782S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	FC-L	All
40107S G4S Linker (E702V) Rabbit mAb (Alexa Fluor® 700 Conjugate)	100 µl	76,000	FC-L	All
17621S G4S Linker (E702V) Rabbit mAb (Biotinylated)	100 µl	76,000	FC-L	All
38907S G4S Linker (E702V) Rabbit mAb (PE Conjugate)	100 µl	76,000	FC-L	All
57710S Whitlow/218 Linker (E3U7Q) Rabbit mAb	100 µl	68,000	FC-L	All
66159S Whitlow/218 Linker (E3U7Q) Rabbit mAb (BSA and Azide Free)	100 µl	134,000	FC-L	All
55809S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 488 Conjugate)	100 µl	76,000	FC-L	All
25186S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 532 Conjugate)	100 µl	76,000	FC-L	All
85651S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 555 Conjugate)	100 µl	76,000	FC-L	All
61465S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 594 Conjugate)	100 µl	76,000	FC-L	All
69310S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 647 Conjugate)	100 µl	76,000	FC-L	All
72160S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Alexa Fluor® 700 Conjugate)	100 µl	76,000	FC-L	All
32523S Whitlow/218 Linker (E3U7Q) Rabbit mAb (Biotinylated)	100 µl	76,000	FC-L	All
62405S Whitlow/218 Linker (E3U7Q) Rabbit mAb (PE Conjugate)	100 µl	76,000	FC-L	All

CST最新情報



AACR 2023の要約

「がん科学と医学の促進」に焦点をあてたAACR 2023では、がん研究と医薬品開発における重要な進歩が注目されました。市場初のCARリンカー抗体や、Agilent Technologies社と行った弊社の発表内容、学会で得た重要な知見などをまとめたブログをご覧ください。

ブログはこちら：

www.cst-science.com/23npi-jul-info1




ストレス顆粒：神経変性疾患におけるその機能と役割

ストレス顆粒 (SG) は、ストレスに対する重要な細胞の防御機構です。SGがどのように形成され、どのタンパク質がSGに局在しているかについてや、様々な神経変性疾患における役割を紹介します。

ブログはこちら：

www.cst-science.com/23npi-jul-info2



がん療法の開発のためのソリューション

ワークフローや、プラットフォームとテクノロジー、がん研究の流れに沿ったソリューションで医薬品開発を加速してください。リソースページでは、効果的な新しい治療法につながる可能性がある標的の情報を紹介しています。

詳細はこちら：

www.cst-science.com/23npi-jul-info3



CSTテクニカルセミナーシリーズ

無料オンラインセミナーを毎月実施しています。弊社の技術者が培った、実験をスムーズに進めるためのコツをお伝えします。どなたでもお気軽にご参加ください。

詳細はこちら：

www.cst-science.com/23npi-jul-info4



種交差性：H-ヒト、M-マウス、R-ラット、Hm-ハムスター、Mk-サル、Mi-ミンク、C-ニワトリ、Pg-ブタ、B-ウシ、Dm-キイロショウジョウバエ、Z-ゼブラフィッシュ、X-アフリカツメガエル、Sc-出芽酵母、Hr-ウマ、Ce-線虫、GP-モルモット、Dg-イヌ、Rab-ウサギ、Vir-ウイルス、All-すべて
括弧付きの動物種は、配列が100%相同であるため反応すると推定されます。種交差性は、ウェスタンブロッティングで確認しています。
アプリケーション：WB-ウェスタンブロッティング、IP-免疫沈降、IF-IC-免疫蛍光細胞染色、IF-P-免疫蛍光染色 (パラフィン切片)、IF-F-免疫蛍光染色 (凍結切片)、IHC-P-免疫組織化学染色 (パラフィン切片)、IHC-LB-Leica® Bond™使用免疫組織化学染色 (パラフィン切片)、IHC-F-免疫組織化学染色 (凍結切片)、FC-FP-フローサイトメトリー (固定、透過化)、FC-L-フローサイトメトリー (生細胞)、ChIP-クロマチン免疫沈降、ChIP-seq-クロマチン免疫沈降シーケンス、C&R-CUT&RUN、ELISA-ELISA、E-P-Peptide ELISA、R Dot Blot-RNAドットプロット、FUNC-Functional

本パンフレット中の情報は2023年6月現在のもので、予告なく変更される場合がありますので予めご了承ください。

© 2023 Cell Signaling Technology, Inc. Cell Signaling Technology, CST, PathScan, and XP are trademarks of Cell Signaling Technology, Inc. Alexa Fluor® is a registered trademark of Life Technologies Corporation. All other trademarks are the property of their respective owners.

For Research Use Only. Not for Use in Diagnostic Procedures. 23-PRG-26351-FLY-JP

輸入販売元



〒101-0047

東京都千代田区内神田1-6-10 笠原ビルディング10階

www.cellsignal.jp

取扱店